


EUV microscopy - a user's perspective

Dimitri Scholz

EUV, 09.11.2011



Begin *ab ovo*

- Simple approaches direct to the goal
- Standard methods of sample preparation

Brainstorming before Experiments

- One hour thinking saves months of frustrations

Light-Electron marriage

Glorious history of microscopy


**Leeuwenhoek
Microscope
(circa late 1600s)**



Cuff-style microscopes,
popular in the mid-1700s,
were the first to provide
ease of use and accurate
focusing mechanisms.



**Pacini
Compound
Microscope
(circa 1845)**



**Leitz
Photomicrographic
Apparatus
(circa 1910)**



Nikon 80i (2008)



Make me a confocal picture!

Electron
microscopy

Light microscopy

Transmission light

Fluorescent light

Structured
Illumination

Bright field

Polarized light

Nomarski contrast (DIC)

Epi-fluorescent

Confocal

TIRF

SPIM

Two-Photon

Decon-
volution

Single pinhole

Parallel array

Spinning disc

Total
Internal
Reflection
Fluorescence

Single
Plane
Illumination
Microscopy

Electron microscopy

Light microscopy

TEM

Biomed

Materials

W or LaB6
cathode

FEG
cathode

80 kEV

200-300
kEV

**Resolution
1 nm**

**Resolution
0.1 nm**

SEM

Table - top

W cathode

5-30 kEV

**Resolution
ca. 5-20 nm**

High spec

FEG cathode

80-300 kEV

Resolution ca. 1 nm

BSE

E₂

E₂

FIB

X-ray

Fluor

BSE

EDAX

Frozen

Environ
mental

Low vac


Please no Eierliegende Wollmilchsau!

Transmission light microscopy

Mean diameter of mouse capillaries

Control $2.6 \pm 0.07 \mu\text{m}$

Epo +/- $3.2 \pm 0.06 \mu\text{m}$






Transmission light microscopy


Growth of collateral arteries continues under zero shear stress

3d ischemia+14d reperfusion

3d ischemia



Polarized light, DIC microscopy



Transmission light microscopy, DIC



Fe/Au Nanoparticles ca. 1 μm

Polarised Light



Condenser Oil NA 1.4

Objective 100x TIRF NA 1.49


Best possible resolution 2011



Fluorescence \neq Confocal



Microtubules: Imunolabeling against α -Tubulin



10μm



3D-Reconstruction of ca. 100 confocal optical Sections

Zebrafish retina

FLUOVIEW® FV10i



Low spec 2009 > High spec 1999



Image Resolution

Resolution: the minimum separation (s) necessary between two point objects in the sample so that they can be distinguished as separate





Image Resolution is limited by two factors:

1. Optical Resolution
2. Detector Resolution

Spatial (XY) resolution



Spatial resolution = the minimal distance between two objects to recognize them as separated



*Founder of
modern
microscopy*

$$R = \frac{1.22 \lambda}{NA_{\text{obj}} + NA_{\text{cond}}} \approx \frac{0.61 \lambda}{NA} \approx \frac{\lambda}{2}$$

Diffractional resolution limit: point objects are detected as point spread functions (PSF)

R

Calculated XY resolution for popular objectives






$$R = \frac{0.61 \lambda}{NA}$$

Lens	NA	400 nm	500 nm	600 nm	700nm
10	0.3	813	1017	1220	1423
10	0.40	610	763	915	1068
20	0.50	488	610	732	854
20	0.75	325	407	488	569
40	0.60	407	508	610	712
40	0.75	325	407	488	569
40	0.95	257	321	385	449
40	1.00	244	305	366	427
40	1.30	188	235	282	328
63	1.20	203	254	305	356
63	1.40	174	218	261	305
100	1.40	174	218	261	305
100	1.49	164	205	246	287

Resolution limitation by detector

Nyquist–Shannon sampling theorem:

Converting from an analog signal (sound or image) to digital, the sampling frequency must be greater than twice the highest frequency of the input signal in order to be able to reconstruct the original perfectly from the sampled version.



3:1

Agreed:
2.4:1

Not number of pixels, but their size, defines resolution

Lens	NA	500 nm	2x	Digital camera pixel size				
				3.4 μm	6.45 μm	8 μm	13 μm	16 μm
10	0.30	1017	508	340	645	800	1300	1600
10	0.40	763	381	340	645	800	1300	1600
20	0.50	610	305	170	323	400	650	800
20	0.75	407	203	170	323	400	650	800
40	0.60	508	254	85	161	200	325	400
40	0.75	407	203	85	161	200	325	400
40	0.95	321	161	85	161	200	325	400
40	1.00	305	153	85	161	200	325	400
40	1.30	235	117	85	161	200	325	400
63	1.20	254	127	54	102	127	206	254
63	1.40	218	109	54	102	127	206	254
100	1.40	218	109	34	65	80	130	160
100	1.49	205	102	34	65	80	130	160


Paradox: High resolution detection is more important for low power objectives.

Rat ventricular myocyte labeled for vinculin (green) and caveolin-3 (red), with colocalized voxels white

A

voxel size 100 x 100 x 250 nm,
which satisfies the Nyquist criteria


The colocalization 19%




B

voxel size 400 x 400 x 400 nm
under-sampled image


The colocalization 61%





C






D




Metric scale and the ultrastructure



From Molecular Biology of Cell. 4th or 5th edition






10 μm

5 x 15 details
for a platelet

Less than 2 for
Z-dimension!



1 μm

Light microscopy vs.

Electron microscopy

Advantages:

- 1) Simple
- 2) Live cells
- 3) Large area
- 4) Multiple labeling

Disadvantages:

- 1) Diffraction-limited XY-resolution ca.
200 nm \approx 50 lines for 10 μ m cell
- 2) Poor Z-resolution $>$ 600nm mostly top vs. bottom
- 3) Poor recognition of organelles

Advantages:

- 1) High resolution
- 2) Structural recognition


Disadvantages:

- 1) Vacuum kills: no live cells
- 2) Only thin sections or surface
= limited info
- 3) Poor immunolabeling
- 4) Multiple labeling difficult
- 5) Expensive



Super-resolution or
correlative microscopy?

STED = Stimulated Emission Depletion



Original fluorescence
excitation spot,
Depletion Laser off

200 nm




Two superimposed beams:


- Excitation laser Pulsed ($<10\text{ps}$) $\rightarrow 635\text{nm}$
- Depletion laser: Pulsed donut-shaped red-shifted (IR, $200\text{-}300\text{ps}$)
- Dyes that can be depleted effectively with low re-excitation \rightarrow ATTO 647N or ATTO655




STED spot size reduction principle




Depletion Power: 1




STED spot size reduction principle




Depletion Power: 2




STED spot size reduction principle




Depletion Power: 3




STED spot size reduction principle




Depletion Power: 4




STED spot size reduction principle




Depletion Power: 5




STED spot size reduction principle




Depletion Power: 6




STED spot size reduction principle



Depletion Power: 7




STED spot size reduction principle




Depletion Power: 8

XY Resolution
in STED is mainly
determined by depletion
power




Myosin in Mausmuskelfaser

confocal



STED



Superresolution: comparison

SIM





Structured Illumination Microscopy

STED


Stimulated Emission Depletion

PAL-M

PhotoActivated Localization Microscopy




Confocal Resolution



10 μm


If 40 nm instead of 220,

Ca. 30 x 80 details for a
platelet



1 μm

CLEM = Correlative Light Electron Microscopy



**Multicolor and Electron Microscopic Imaging of
Connexin Trafficking**
Guido Gaietta, *et al.*
Science **296**, 503 (2002);
DOI: 10.1126/science.1068793

C L E M = Correlative Light Electron Microscopy

Transfection Cx43-tetracysteine

Live cells: ReAsH (A)


Fluorescent microscopy

Fixation

DAB+ UV O^*

Epon

EM (B, C)



CLEM = Correlative Light Electron Microscopy

Fixierung

Anti- β -tubulin



Secondary*quantum dot_655

DAPI

Fluoreszenzmikroskopie (A)


Fixierung - Epon

EM (B)



Up Fluorescence, down SEM

JEOL JASM-6200 Scanning Electron Microscope ClairScope 2010



We develop an X-ray-Fluorescence combination

	Practically achievable XY resolution, nm	Features, Limitations
Fluorescence microscopy, including confocal and multi-photon	200 nm	<ul style="list-style-type: none"> - simple, suitable for live cells - multiple labeling, large field of view - suitable for fast acquisition (30 fps and more) - low resolution
Electron microscopy	1 nm (bio-EM) 10 nm (immuno-EM)	<ul style="list-style-type: none"> - vacuum: unsuitable for most live cell studies - thin samples only (70-300nm): reduced information - expensive sample preparation - immunolabeling difficult and decreases the resolution
Super-resolution fluorescence: STORM, PALM, STED	40 nm	<ul style="list-style-type: none"> - slow: unsuitable for many live cell studies - require photo switchable fluorescent dyes or proteins
Near UV microscopy	100 nm	<ul style="list-style-type: none"> - requires special expensive lenses - poor signal/noise ratio - only 2-fold gain of resolution
EUV microscopy now	40 nm	<ul style="list-style-type: none"> - requires special light sources and X-ray optics - narrow field of view