

Compact EUV Source for Metrology and Inspection

Klaus Bergmann, Jochen Vieker, Alexander von Wezyk
2015 EUV Source Workshop, 10.11.2015, Dublin

Overview

- Introduction
- Xenon based EUV Source FS5420
- Consideration on 6.x nm emission
- Recent progress on lifetime and power scaling
- Dose control of brilliance

EUV Source Activities at Fraunhofer ILT





AIxUV
EUV-Technology

ri research
instruments



PHILIPS

Tin vacuum arc



PHILIPS

Sources for Aachen University

- Water window microscopy (2.9 nm)
- Nano patterning
- Coherent diffraction imaging
- EUV microscopy
- Photo electron spectroscopy



TOS
RWTH
RHEINISCH-
WESTFÄLISCHE
TECHNISCHE
HOCHSCHULE
AACHEN
LEHRSTUHL
FÜR TECHNOLOGIE
OPTISCHER SYSTEME



EUV
RWTH
RHEINISCH-
WESTFÄLISCHE
TECHNISCHE
HOCHSCHULE
AACHEN
CHAIR FOR
EXPERIMENTAL PHYSICS
OF EXTREME ULTRAVIOLET

XTREME
technologies

USHIO
GROUP

v:\a11.2015.tagung.euv\2015.dublin>S32_dublin2015.pptx

Hollow Cathode triggered Pinch Plasma



Fraunhofer ILT Source : FS5420

■ Standard Mode

Inband Power : $> 20 \text{ W}/2\pi\text{sr}$

EUV pulse energy : 2,2 mJ/sr

typ. repetition rate : 1500 Hz

ave. peak brightness : $8 \text{ W}/\text{mm}^2\text{sr}$

■ High Pulse Energy Mode

Inband Power : $< 10 \text{ W}/2\pi\text{sr}$

EUV pulse energy : $> 4,0 \text{ mJ/sr}$


typ. repetition rate : $< 400 \text{ Hz}$



EUV source including power supply, control unit & chillers

HMI - Software

- repetition rate and pulse energy selectable in pre-defined ranges
- gas flow automatically adjusted according breakdown delay
- option for soft ramping up and down of input power
- display of measured frequency and power
- measuring and display of missing pulses (Quality Factor)
- option for EUV power monitoring and dose-control
- fully accessible via Ethernet / Fieldbus




v:\a11.2015.tagung.euvl\2015.dublin>S32_dublin2015.pptx


Emission Characteristics

- EUV source dimensions:
 - diameter (FWHM) : $< 300 \mu\text{m}$
 - length (FWHM) : $\sim 4 \text{ mm}$
- stability of diameter : $\sim 4 \%$
- spatial stability : $< 7 \mu\text{m}$
- pulse duration : $50 - 100 \text{ ns}$
- power stability : $1-3 \%$
(100 pulses moving average
over 10 minutes without dose control)

Typical F55420 emission spectrum






Axial and radial emission profile



Emission around 6.x nm


- Contributions from 4d-4f transitions around KrX
- Emission between 6-7 nm at 4,4 kW input power : $\sim 15\text{W}/2\pi\text{sr}$
- Suitable for irradiation damage studies or optics characterization




v:<a11.2015.tagung.euvl2015.dublin>S32_dublin2015.pptx

Power Scaling: FS5420 in 40 W - Operation



- average input power:
13 kW
- repetition rate:
2500 Hz
- peak brightness:
 $\sim 12 \text{ W/mm}^2\text{sr}$
- pulse-to-pulse standard deviation:
6,8 %



Power Scaling: FS5420 in 55 W - Operation



Source Development Aspects




- better understanding of origin of decrease of CE and performance during operation
- reduction of erosion by choosing sputter resistive material
- current achievement for the cathode: **>800 Mshot** and still running




Comparison of Mo and W-based cathode after same number of Mshots

Advanced Triggering



Conventional triggering

Advanced triggering




- More sophisticated $U_{\text{trigger}}(t)$ at trigger electrode
- Use of additional potential at intermediate electrode $U_{\text{IP}}(t)$

v:<a11.2015.tagung.eu>v12015.dublin>S32_dublin2015.pptx

Increase of operation window

- larger process window for Xe- flow due to **advanced triggering**
- higher tolerance towards electrode erosion
(tested electrodes had >150 Mshot)
- higher CE (~0.7 %) at lower gas flows achievable
- 13.5 nm inband power at 6 kW input:
 $\sim 40 \text{ W}/(2\pi\text{sr})$
- 13.5 nm inband brightness at 6 kW input:
 $\sim 12 \text{ W}/(\text{mm}^2\text{sr})$




Former process window

Increase of maintenance interval

- wear of IP (Intermediate Plate) is mainly determining the source performance
- Xe flow increases with wear of IP
- higher flow lead to lower CE and reduces the working range
- advanced triggering allows operation at lower Xe flow with increase of CE
- **1 Gshot** maintenance interval with advanced trigger expected



Intermediate Plate (IP)



Demonstration of Dose Control of Inband-Power

Dose control at Philips-EUV Xenon source for ASML alpha-tool




Figure 4: Effect of the dose control system. (a) Dose values over 60000 subsequent pulses. (b) Histograms over about 50000 shots with and without dose control system. The dose stability is improved by a factor of 10.

Ref. : Pankert et al., SPIE 2005


- measurement of inband EUV for each pulse
- assumption of linear behavior of EUV pulse energy with input pulse energy
- variation of input pulse energy for each pulse
- expected dose stability within window of N pulses (σ : standard deviation)

without dose control
with dose control

: σ / \sqrt{N}
: σ / N


Dose Control of Brilliance : Task

- Potential improvement of $1/N^{1/2}$ allows for relaxed conditions for source stability σ_B
- Look for suitable parameter with high correlation, ρ , with source brilliance
- Built up of pulse-to-pulse metrology for this parameter (e.g. total inband photon flux)



Dose Control of Brilliance : Requirement on Correlation


without dose control




- Simulation for Gaussian Distributions $\langle x \rangle$, $\langle y \rangle$ with $\sigma = 10\%$ each and average values of 1.0
- x measured parameter, e.g. inband photon flux
- y controlled parameter, i.e., source brilliance
- Integration over $N = 100$ pulses
- improved dose stability requires $\rho > 0.5$

with dose control


$\rho = 0.3$



$\rho = 0.7$





$\rho = 1.0$



v:\a11.2015.tagung.euvl2015.dublin>S32_dublin2015.pptx


Dose Control of Brilliance : Correlation Parameter

- Proof of concept experiment at 1 kHz & 4 J input power
Fast EUV-inband Cameras at 2,5 Hz pulse picking
- Metrology:
 - Camera on-axis
 - Camera off-axis
 - EUV Energy Monitor off axis
- Evaluation at fixed pixel-position
- High sensitivity setup without Zr-Filter at off-axis camera
-> some visible light may still be interfering
- High Correlation Coefficient $\rho > 0,85$ for Spot with 60 μm diameter



Single-shot Camera images of pinch:
On-Axis

Off-Axis



Linear Brightness – Pulse energy dependence allows
active control of dose collected from small NA

Summary

- New generation FS5420 EUV source is available at Fraunhofer ILT
- Improvements in operation by software, e.g. error detection and automatic adjustment of operation parameters
- Ongoing basics investigations on lifetime and power scaling
- Major improvements have been demonstrated in key experiments by means of advanced triggering
- Next generation **FS5440** – 40 W EUV source in development:
Combining long-life electrodes with high CE of new trigger system