

# Laboratory-based Photoemission Spectro- microscopy at 71.7 eV and 113.7 eV for Studies of Complex Materials

*Daniel Wilson, Christoph Schmitz, Denis Rudolf, Sally Rieß, Martin  
Schuck, Carsten Wiemann, Astrid Besmehn, Hilde Hardtdegen,  
Detlev Grützmacher, Claus M. Schneider, F. Stefan Tautz,  
and Larissa Juschkin*

**8 December, 2015**

## Motivation: Intense EUV Light Sources

# New Journal of Physics

The open-access journal for physics

## Photoemission spectroscopy—from early days to recent applications

**Friedrich Reinert<sup>1,3</sup> and Stefan Hübner<sup>2</sup>**

“Furthermore, one can hope that in the future there will be intense and continuous light sources available for the laboratory, which will end the necessity to move to synchrotron radiation facilities for numerous photoemission applications.”

Friedrich Reinert, 2005, doi:10.1088/1367-2630/7/1/097

# Motivation: Intense EUV Light Sources



- Standard sources for UPS limited to 40.8 eV, but problems arise:
  - e.g. final state effects
- For XPS, X-ray tubes are used
  - Long integration times (small interaction cross sections)
  - No UPS

	Laboratory Light Sources	Photon Energy Range
UPS	UV-Laser	4.66 eV
	Hg-Lamp	4.9 eV
	He I/II	21.2 eV / 40.8 eV
UPS & XPS	HHG	~ 100 eV
	DPP-Source	Up to 620 eV
XPS	X-ray Tube	1488 eV ( $K_{\alpha}$ )

# Motivation: Intense EUV Light Sources

- Standard sources for UPS limited to 40.8 eV, but problems arise:
  - e.g. final state effects
- For XPS, X-ray tubes are used
  - Long integration times (small interaction cross sections)
  - No UPS
- HHG sources suffer from extreme space charge effects
  - Short EUV pulses -> extremely high PE densities
  - Coulomb repulsion at crossings
  - Temporal, lateral & energy broadening

UPS & XPS	HHG	~ 100 eV
	DPP-Source	Up to 620 eV





Time-of-flight-photoelectron emission microscopy on plasmonic structures using attosecond extreme ultraviolet pulses, Chew et al. 2012

Integration time: 1.5 h,  $E_{\gamma} = 93$  eV


# Setup

- $E_\gamma = 71.7 \text{ eV}$
- $E_{PP} = 1.32 \frac{\text{mJ}}{2\pi \text{ sr}}$
- $\bar{I} = 2.64 \frac{\text{W}}{2\pi \text{ sr}}$
- On sample:
  - $1.8 \cdot 10^{11} \text{ photons}/(\text{mm}^2 \text{ s})$



# The NanoESCA

- Analysis of chemical and electronic states in **complex** and **nanostructured** material
- Electrostatic **PEEM** lenses and a double-pass hemispherical analyser
- “Resolutions in **imaging ESCA** of 650 nm in the laboratory and **150 nm** at the **synchrotron** are simply unique.”  
<https://www.elettra.trieste.it/elettra-beamlines/nanoesca.html>
- NanoESCA used here is identically constructed






# Modes of the NanoESCA

Where are the interesting features?

Which elements are there and how much?


What is the elemental map?

**PEEM Mode**




**Photoelectron Emission Microscope mode**

**Spectroscopy Mode**



**Small Spot Spectroscopy mode**


**Imaging ESCA Mode**




**Imaging ESCA-Mode**

# EUV spectrum (oxygen) behind Bragg mirrors

emission spectrum of oxygen plasma



- Bragg mirrors from **optiXfab.**





- Doppler broadening<sup>1</sup>: ~ 13 meV
- $O^{5+}$  separation: 60 meV
- Distance to  $O^{4+}$ : < 400 meV
  - Effect from satellite does not affect measurements

<sup>1</sup>Saha LTE spectrum from NIST database:  
 $T_e = 10 \text{ eV}$ ,  $n_e = 10^{18} \text{ cm}^{-3}$



## PEEM mode: Features of GeSbTe (GST)




# Spectroscopic opportunities in the lab

# Study of epitaxially grown crystalline hexagonal Ge-Sb-Te (GST<sub>1:2:4</sub>) islands on Si (111).


## Focus: Surface oxidation

- a) EF-PEEM image at 5.2 eV (Hg)
- b) EF-PEEM images from 39 – 41 eV (Te 4d) show regions of high Te concentration
  - $E_{\gamma} = 71.7 \text{ eV}$
  - $t_{integration} = 1 \text{ s per image}$
- c) Spectrum extracted from images
  - Black: Oxidized Sb peak visible
  - Blue: sputtered/annealed sample
  - Red: Standard XPS with Al-K $\alpha$
- Standard XPS lead to same results




# Resolution in energy filtered (EF) PEEM mode

- Sample: chessboard calibration
- Au fields  $1 \times 1 \mu\text{m}^2$ 
  - a) EF-PEEM:  $E_{\gamma} = 71.7 \text{ eV}$
  - b) EF-PEEM:  $E_{\gamma} = 5.2 \text{ eV}$
  - c) Line scans for each picture
- Contrast is low for  $E_{\gamma} = 71.7 \text{ eV}$ 
  - Core levels of Si & Au out of range
- Resolution is still comparable:
  - $\Delta x = 320 \text{ nm}$  (HG)
  - $\Delta x = 424 \text{ nm}$  (EUV 71.7 eV)



## k-Space measurements with 71.7 eV



- Maximum field of view:  $5.6 \text{ \AA}^{-1}$ 
  - 2.4 Brillouin zones of Au (111)
  - Six fold symmetry
- Energy resolution  $\sim 400 \text{ meV}$  from electron analyzer
  - Two emission lines at  $17.3 \text{ nm}$ 
    - approx.  $15 \text{ meV}$  linewidth
    - $60 \text{ meV}$  distance






# k-Space measurements with 71.7 eV

- Maximum field of view:  $5.6 \text{ \AA}^{-1}$ 
  - 2.4 Brillouin zones of Au (111)
  - Six fold symmetry
- Energy resolution  $\sim 400 \text{ meV}$  from electron analyzer
  - Two emission lines at  $17.3 \text{ nm}$ 
    - approx.  $15 \text{ meV}$  linewidth
    - $60 \text{ meV}$  distance
- Surface states occur only at the sharp transition from solid material that ends with a surface



# High flux from Xe 10.9 nm UTA


- Adding Argon to Xenon results in enhanced emission at 10.9 nm (113.7 eV) by suppressing lower and higher wavelengths.
- 2x small-bandwidth ML-mirrors
  - $\Delta\lambda_{FWHM} = 0.1 \text{ nm}$
  - $\Delta E_{FWHM} = 1 \text{ eV}$
  - $R_{peak} = 13 \%$











# Debris analysis: EUV source exposure sample

- Sample has been exposed to a LPP EUV source
- PES ( $E_{\gamma} = 113.7$  eV) shows one broad peak at  $E_{bin} = 75$  eV
- Elements found fit perfectly to EUV source & sample history
  - Cu, Al
- Measurement time (complete range): 10 min



# Summary & Outlook

- High power EUV gas discharge sources are well suited for UPS/XPS/PEEM
  - High flux
  - Extreme narrow emission lines
  - Reasonable pulse duration
  - Tunable spectrum
  - Highly surface sensitive
  - Multiple beamlines possible
- Future outlook:
  - More sophisticated optics
  - Improvement of EUV emission characteristics



## Acknowledgements

- Christoph Schmitz
- Denis Rudolf
- Carsten Wiemann
- Prof. L. Juschkin
- Prof. C. M. Schneider
- Prof. F. S. Tautz
- GST Samples
  - Sally Rieß
  - Martin Schuck

