


Exploring Radiance

Fergal O'Reilly, Tony Donnelly, Mateusz
Olszewski, Isaac Tobin, Pauline
O'Callaghan, Tom McCormack,
Gladson Joseph

Soft X-ray Nanoscopy




Soft X-ray Imaging of particular interest if **SPATIAL RESOLUTION better than 50 nm**

Larabell et al: Biol. Cell (2017) **109**, 24–38 DOI: 10.1111/boc.201600044

Imaging Challenge

CCD Cameras are Soft
X-ray imaging
Gold Standard



Smallest **Pixels** > **13 μm**

For 50 nm resolution **object pixels** < **15 nm**
(2D Nyquist Resolution)

So **magnification** > **1000**


Affordable (< €50k) CCD has 25 mm X 25 mm size

At > x1000 mag that means **FOV < 25 μm** field of view

Rayleigh Resolution

For 50 nm resolution at 2.5 nm – **NA > 0.03** (> **0.003 steradians**)

$$\text{Resolution} \sim 1.22\lambda / (\text{NA}_{\text{illumination}} + \text{NA}_{\text{objective}})$$



$$\text{Coherence} = \text{NA}_{\text{illumination}} / \text{NA}_{\text{objective}}$$

Photon Flux


- At Sample:
 - Ideally want >1000 photons/pixel/image
 - This is what the synchrotron systems can achieve
 - Fills the CCD pixel well
 - Gives highest possible contrast

Average Radiance at Sample

- 2D Images taken over second timescales
- $$\text{Radiance} = \frac{\text{SourceFlux} \times \text{SystemTransmission}}{(\text{FOV Diameter}^2 \times \text{Objective Solid Angle})}$$
- *Peak radiance is only of interest if you are looking at fast phenomena*

Source Size and Collection

Source Diameter Reduced by 2




Collect 4x photons from source

Laser Plasmas & Soft X-rays

$$T_e(\text{eV}) \approx (\lambda^2 \text{ PowerDensity})^{3/5}$$

(Columbant & Tonan, JAP 44 1973)




$$\text{Average charge} \approx 0.67 (AT_e)^{1/3}$$

$$n_e \sim 10^{19} - 10^{21} \text{ cm}^{-3} (n_{ec} \sim 10^{21} / \lambda^2 \text{ cm}^{-3})$$

$$\text{Expansion velocity} \approx 10 - 100 \text{ } \mu\text{m/ns} \\ (\approx 10^6 - 10^7 \text{ cm/s})$$

Soft X-ray/EUV Spectrum



For an optically thin plasma:


$$P_{lines} : P_{recomb} : P_{brem} = 100 : 10 : 1$$

Plasma Size Experiments


- Requirements:
 - to **image plasmas** $< 5\mu\text{m}$
 - simultaneous Soft X-Ray spectroscopic study
- Objective:
 - Explore the size and shape of the Soft X-Ray plasmas with different target materials
 - Understand the Soft X-Ray emission volume from different materials at different conditions
 - **Optimise Radiance** at particular wavelengths

Current State of the Art

- Hans Hertz $\sim 14\mu\text{m} \times 20\mu\text{m}$ – measured with condenser zone plate and multilayer optic



Martz et al November 1, 2012 / Vol. 37,
No. 21 / OPTICS LETTERS



UCD Expertise

- Plasma facing optics
- Can use very high laser illumination NA
- So potentially very small plasmas
- Of potentially any target
- Can use inexpensive industrial lasers


High Resolution Soft X-ray Source Imaging



Facilities at UCD SPEC Lab




Lasers



Main chamber and Spectrometers

800 mJ, 7 ns, 1064 nm M2 ~3.5
500 mJ, 170 ps, 1064 nm M2 ~5
50 mJ, 20 ps, 1064 nm M2 ~5


Previous Imaging



Pinholes have been used before

Miyazaki et al, <http://www.euvlitho.com/2015/S53.pdf>





Experiment




Convolve Gaussian Source and Airy

$$\begin{aligned}
 & \text{Gaussian Source} * \text{Airy} = \text{Image FWHM} \\
 & (\text{Magnification} * \text{Source FWHM})^2 + \text{Airy FWHM}^2 = \text{Image FWHM}^2
 \end{aligned}$$


Filtered Spectrum




Experiment Setup



Pinhole at lower magnification



Pinhole at higher magnification



Spectroscopy setup

Pinhole Comparison with Mo plasma




Laser Energy 800mj

200mj




50mj

Pinhole
Diameter




10 μ m




5 μ m



2 μ m




Energy Scan



Images are Energy Normalised

Silicon



Magnification = 94.

Filter = Si₃N₄ and Ag

Plasma Side View - Mo

FWHM $\sim 5 \mu\text{m}$

Image Number:1783






Image Number:1783



FWHM $\sim 16 \mu\text{m}$

Image Number:1734








Image Number:1734




Plasma Front View - Mo




Burning Through




Less Energy = Smaller Plasma



Conversion Efficiency




Smaller Not Always Better



Conclusion

- We have demonstrated high CE from sub 10 μm plasmas
- Much more space to optimise into



Acknowledgments



Erasmus Mundus
A programme of the European Union

We would like to acknowledge the kind support of

All in the UCD Spec Group,
and the Mechanical and Electronic Workshops, and the
Administrative Staff of UCD School of Physics

EXTATIC - Erasmus Mundus Programme

Thank you for listening