

Broad-band spectral characterization of plasma based sources of soft X-ray and VUV radiation



2017 Source workshop, Nov. 6-8, Dublin, Ireland

Applications = motivation




- Search for efficient plasma-based sources of VUV radiation
- High-harmonic generation
- In-band vs out-of-band emission measurements for EUVL sources

1. BROADBAND SPECTROSCOPY

AGS = amplitude grating spectrometer

Important features:


- ★ No focusing elements
- ★ Flat field
- ★ Robust
- ★ Relaxed grating specs
- ★ Very compact



Advantages of the grazing incidence (GI) operation

1

Effective period reduction: $d_{\text{eff}} = d \sin \theta$




~100-1000 l/mm gratings can be used

2

GI illumination with a confined beam -> large number of grating periods in the diffraction, increases resolution

Diffraction orders



$$b = fp$$


Deep grooves - light reflected only from top surface

Diffraction efficiencies predicted by scalar diffraction:

$$I_m = I_{inc} D_m^2 R_m$$

$$D_m = f \operatorname{sinc}(\pi m f)$$

for details see: [JOSA 45\(9\), 756 \(1955\)](#)



even orders suppressed when $f = 0.5$

Handbook diffraction principles


Main gratin equation:

$$\sin \psi + \sin \varphi = m\lambda/p$$

Angular dispersion:

$$D = d\varphi/d\lambda = m/(p \cdot \cos \varphi)$$

D significantly increases at $\varphi \rightarrow 90^\circ$



Implementation



View of device with Hamamatsu CCD detector

SPECIFICATIONS:

Grating period	3 μm
Grazing angle	6 deg
Slit width	80 μm
Spectral region	6-200 nm
Resolution, $\lambda/\delta\lambda$	up to 50
Detector	CCD*

*CCD options

line Toshiba 1304AP, 3600 8 μm pixels

array Hamamatsu S7030/S7031, 1024x122 24 μm pixels

array ANDOR DX440-BN, 2048x512 13.5 μm pixels



Spectral resolution

Calculated for the given hardware configuration (see prev. slide) + using typical values:

Source - spectrometer distance 0.5 - 1 m

Source size 0.1 - 1 mm

reminder for order sign notation




Wavelength calibration

Radiation source: laser-produced plasma with light (low-Z) target materials (LiD, Be, C,...)


Registration system: Back-side illuminated CCD camera ANDOR DX440-BN

NIST tabulated data for optical transitions of light elements is used for the wavelength calibration


LiD



Be



C




Wavelength calibration

Radiation source: laser-produced plasma with light (low-Z) target materials (LiD, Be, C,...)

Registration system: Back-side illuminated CCD camera ANDOR DX440-BN



NIST tabulated data for optical transitions of light elements is used for the wavelength calibration



High resolution and moderate resolution

Comparison between +1 and -1 orders of diffraction

Radiation source: laser-produced plasma with LiF target




evidently, “ $m = -1$ ” lines are better resolved

Contribution of high diffraction orders

Radiation source: laser-produced plasma with Sn target

Registration system: Back-side illuminated CCD camera ANDOR DX440-BN




Major contribution is due to the main Sn peak at 13.5 nm

Contribution of high diffraction orders

Radiation source: laser-produced plasma with Sn target

Registration system: Back-side illuminated CCD camera ANDOR DX440-BN


+Spectral filter: free-standing Zr/Si multilayer with bandpass centered at 13.5 nm



⇒ Diffraction efficiencies can be estimated for further spectral postprocessing

2. HIGH-RESOLUTION SPECTROSCOPY

Rowland and off-Rowland: how it usually works




R1 = Rowland - detector has to be bent
-> not compatible with standard CCD

R2 = off-Rowland - detector is out of the focal plane
-> high resolution achieved only at around the focal spot

What if we shift the entrance slit?

Transformation of the focal curve vs slit position



G - grating; D - detector; RC - Rowland circle




It is possible to design optical scheme with nearly flat field perpendicular illumination of the detector

Old idea revision - see W. Baily, [Philosophical Mag. Ser. 5](#) **15**, 183 (1883).

Implementation for grazing incidence

Test experimental configuration:


Line density	600, 1200 l/mm
Grating radius of curvature	1 m
Working size of the grating	6 mm
Entrance slit width	10 μm
Angle of incidence	4 deg
Slit-grating distance	33 mm
Detector-grating distance	320 mm
Diffraction order	1st



Experimental results

Radiation source:
discharge in argon


Registration system:
Back-side illuminated
CCD camera ANDOR
DX440-BN



Experimental results

Radiation source:
discharge in argon

Registration system:
Back-side illuminated
CCD camera ANDOR
DX440-BN



Summary

- Developed compact grazing-incidence spectrometer for broadband SXR-VUV spectral measurements
 - Spectral range 5-200 nm
 - Spectral resolution up to 50 in the range in 1st diffraction order
 - Spectral resolution up to 100 in the range in -1st diffraction order (but limited spectral range)
- Developed grazing-incidence spectrometer for high resolution measurements in EUV range
 - Spectral range 10-40 nm
 - Spectral resolution up to 600
- For more metrology solutions see www.isteq.nl

Contributing authors



P. S. Antsiferov, A. P. Shevelko, D. B. Abramenko, L. A. Dorokhin, P. V. Krainov,
V. M. Krivtsun, Yu. V. Sidelnikov, V. V. Medvedev, K. N. Koshelev

