


# Compact, bright , Plasma-based EUV Lasers for Metrology

**Jorge J. Rocca**, B. Reagan, Y. Wang,  
D. Alessi, B. Luther, K. Wernsing, L. Yin,  
M. Curtis, M. Berrill, D. Martz, V. Shlyaptsev,  
S. Wang, F. Furch, M. Woolstron, D. Patel,  
M.C. Marconi, C.S. Menoni



**Engineering Research Center for  
Extreme Ultraviolet Science & Technology  
Colorado State University**

**Work Supported by the NSF Engineering Research Centers  
Program and the US Department of Energy**


**Microscopy**



**Nano-printing**




**Analytic nanoprobe**



Applications require EUV laser pulses with pulse high energy (large number of photons), high average power

## Compact atomic soft x-ray lasers







- High pulse energy ( $\mu\text{J}$ - $\text{mJ}$ )
- High monochromaticity ( $\lambda/\Delta\lambda < 10^{-4}$ )
- High peak spectral brightness

Single shot  
imaging

Nanoscale  
ablation

Analytic  
nanoprobes


Error-free  
nanopatterning





# EUV lasers can be created by electron impact excitation of highly ionized atoms in dense plasmas

**Singly ionized** Ar ion, Kr ion lasers in the visible spectral region




Plasma requirements:  
Te ~ 5 eV  
Ne ~  $1 \cdot 10^{14} \text{ cm}^{-3}$


NexTe increases by  $10^7$ - $10^{10}$

Te ~ 100- 1000 eV  
Ne ~  $1 \cdot 10^{19}$  -  $1 \cdot 10^{21} \text{ cm}^{-3}$

**Highly ionized** (8-35 times) in the EUV/SXR spectral region




Laser created plasma



Discharge created plasma

$$h\nu = \Delta E \propto Z^{2.5}$$



**Ionized 20 times**

## Assessing the wavelength extensibility of optical patterned defect inspection

Bryan M. Barnes\*, Hui Zhou, Mark-Alexander Henn, Martin Y. Sohn, and Richard M. Silver Engineering Physics Division, National Institute of Standards and Technology,  
100 Bureau Drive MS 8212, Gaithersburg, MD, USA 20899-8212


### ABSTRACT

Qualitative comparisons have been made in the literature between the scattering off deep-subwavelength-sized defects and the scattering off spheres in free space to illustrate the challenges of optical defect inspection with decreasing patterning sizes. The intensity scattered by such a sphere (for diameters sized well below the wavelength) is proportional to its diameter to the sixth power, but also scales inversely to the fourth power of the wavelength. This paper addresses through simulation the potential advantages of applying shorter wavelengths for improved patterned defect inspection. Rigorous finite-difference time-domain 3-D electromagnetic modeling of the scattering from patterned defect layouts has been performed at five wavelengths which span the deep ultraviolet (193 nm), the vacuum ultraviolet (157 nm and 122 nm), and the extreme ultraviolet (47 nm and 13 nm). These patterned structures and defects are based upon publicly disclosed geometrical cross-sectional information from recent manufacturing processes, which then have been scaled down to an 8 nm Si linewidth. Simulations are performed under an assumption that these wavelengths have the same source intensity, noise sources, and optical configuration, but wavelength-dependent optical constants are considered, thus yielding a more fundamental comparison of the potential gains from wavelength scaling. To make these results more practical, future work should include simulations with more process stacks and with more materials as well as the incorporation of available source strengths, known microscope configurations, and detector quantum efficiencies. In this study, **a 47 nm wavelength yielded enhancements in the signal-to-noise by a factor of five compared to longer wavelengths and in the differential intensities by as much as three orders-of-magnitude compared to 13 nm, the actinic wavelength for EUV semiconductor manufacturing.**

Metrology, Inspection, and Process Control for Microlithography XXXI, edited by Martha I. Sanchez, Vladimir A. Ukraintsev Proc. of SPIE Vol. 10145, 1014516 · © 2017 SPIE · CCC code: 0277-786X/17/\$18 · doi: 10.1117/12.2262191



A fast current pulse compresses a plasma creating a hot and dense column with aspect ratio  $L/d > 1000$




Sequence of on-axis interferograms showing rapid plasma column compression in an Ar capillary discharge

## Demonstration of a Discharge Pumped Table-Top Soft-X-Ray Laser


J. J. Rocca, V. Shlyaptsev,\* F. G. Tomasel,<sup>†</sup> O. D. Cortázar,<sup>‡</sup> D. Hartshorn, and J. L. A. Chilla

*Electrical Engineering Department, Colorado State University, Fort Collins, Colorado 80523*


(Received 31 May 1994)



# Essentially full spatial coherence is achieved increasing the capillary length







Ne-like Ar Capillary discharge 46.9 nm laser

High average power: up to 3 mW

High pulse energy: 0.1 mJ - 0.8 mJ @4 Hz

Narrow spectral bandwidth:  $\Delta\lambda/\lambda = 3 \times 10^{-5}$

Beam divergence:  $\theta = 4.5$  mrad

# CSU Capillary Discharge Laser Technology is today turn-on key



**2005**

- 10 microjoule /pulse
- 0.15 mW average power
- 1-12 Hz repetition rate
- Pulse duration  $\sim 1.5$  ns
- $\Delta\lambda/\lambda < 1 \times 10^{-4}$

*S. Heinbuch, M. Grisham, D. Martz, J.J. Rocca  
Optics Express, 30,2095, (2005)*



**2015**

- **50** microjoule /pulse
- 0.5 mW average power
- 1-10 Hz repetition rate
- Pulse duration  $\sim 1.5$  ns
- $\Delta\lambda/\lambda < 1 \times 10^{-4}$
- Jitter < 2ns









46.9 nm SXR laser

Microscope  
vacuum chamber

## TRANSMISSION CONFIGURATION

Single shot image of  
50 nm nanotubes

2  $\mu$ m

Courtney  
BrewerFernando  
Brizuela




# EUV microscopes captures images of nanostructures with very high resolution

## TRANSMISSION


NA=0.32 (M~1000)

200 nm half period diatom



5 sec exposure

50 nm carbon nanotube




Single shot,  
1.5 ns  
exposure

*C. Brewer, et al, Opt. Lett. 33, 518 (2008)*


## REFLECTION

NA=0.12 and 0.19 (M~250)


5-20 sec exposures



Partially Processed  
semiconductor chip






100 nm thick GaN nanowire  
between Al contacts





Zr surface showing twin in grain


*F. Brizuela et al, Opt. Express 13 :3983, (2005).*



## Magnetic force microscope tip interaction with stray magnetic field




Magnetic field along z




**Effective Spring Constant**  
 $k_{\text{tip}} + k_{\text{force}}$

$$\omega_{\text{res}}^2 = k/m = 1/m (k - \partial F/\partial z)$$






# Real time nano-scale resolution imaging by Fourier Holography




## Numerical back-propagation




Measurements to validate a model that describes the nonlinear interaction between nanopillars

# Applications in dense plasma diagnostics and photochemistry




## Plasma Interferometry


J. Filevich et al PRL 94, 035005 (2005)



M. Purvis et al. Phys. Rev.E, **76**, (2007); **124**, (2010)



## Single photon ionization mass spectrometry




F. Dong et al. J.Chem.Phys **124**, (2006)


F. Dong et al. J.Am.Chem Soc. **131**, (2009)

# Beyond morphology: composition

Condensed phase micro- and nanostructures




Micro-organisms






# EUV laser ablation time-of-flight (TOF) imaging mass spectrometry




1. Focused laser ablates the material
2. Simultaneously ionizes species in the ablation plume
3. Ions are extracted into the mass spectrometer

# 3D imaging EUV mass spectrometry imaging nanoprobe



**EUV Optics**  
Fresnel zone plate 0.16 NA



## **EUV laser**



- Wavelength: 46.9 nm
- Energy per pulse > 10 μJ
- Repetition rate: 12 Hz
- Pulse duration ~1.5 ns

*S. Heinbuch Optics Express  
vol. 13, 4050 (2005)*

Optics engineered by  
W. Chao and E. Anderson at  
Center of X-Ray Optics,  
Lawrence Berkeley Lab.

*E.H. Anderson, IEEE J. Quantum  
Electronics, vol.42, 27, 2006*

# 80 nm lateral resolution was measured on both inorganic and organic samples





***I. Kuznetsov et al. Nature Comm., 6, 6944,( 2015)***


***T. Green, et al., J. of Analytical Atomic Spectrometry, 32, 1092, (2017)***



# Depth resolution of 20 nm is achieved on organic multilayers




SIMS-TOF  
Evans Analytics




I. Kuznetsov et al, Nature Communications, **6**, 6944, (2015)

# 2D Imaging of uranium micron size particles by EUV TOF MSI reveals isotopic content




**Measured:**  
 $^{235}\text{U}/^{238}\text{U}$  ratio of **0.248 ( $\pm 0.0014$ )**  
**Expected value is**  
**0.2513 ( $\pm 0.0003$ ) (U200 standard)**

Collaboration with Pacific  
Northwest Laboratory




*Error free printing*



*M. Marconi, F. Cerrina, et al. (2009)*


**Proof of principle:  
120 nm resolution**







# Error-Free Printing of Periodic metallic structures

Talbot Mask




Print in HSQ



**Grating etched in Au  
500/600 lines/spaces**



# Scaling to shorter wavelengths



# Simulation showed gain-saturated amplification at 13.2 nm in Ni-like Cd can be achieved with ~ 1 J pump

Pre-pulse  
300 mJ, 120 ps



Heating pulse  
1 J, 6 ps






# EUV lasers excited by rapid heating of plasmas with short laser pulses

## Laser Pumping Geometry




Grazing incidence allows for efficient heating of plasma region with optimum electron density




*R. Keenan et al, Phys. Rev. Lett. 94, 103901 (2005) ; B.M. Luther et al, Opt. Lett. 30, 165 (2005); Transient excitation: P. Nickels, V. Shlyaptsev et al. Phys. Rev.Lett. 78,2 748, (1997)*

# High repetition rate table-top EUV lasers in transitions of Ni-like ions down to 10.9 nm




Gain saturated  
operation  
demonstrated



Y. Wang et al, *Phys. Rev. A* **72**, 053807 (2005)



\*D. Martz et al. *Optics Lett.* **35**, 1632 (2010)





- Uses output from 13.2 laser  
*J.J. Rocca et al., Opt. Lett. Vol. 30, 2581 (2005); Y. Wang et al., Phys. Rev. A 72, 5 (2005)*
- Captures images with large field of view with exposures of 20-90 seconds

## TRANSMISSION


- Spatial resolution: better than 38 nm



*G. Vaschenko, et al Optics Letters, vol 31, 1214 (2006)*




LER/CD=0.08



## REFLECTION

- Zone plate microscope for at wavelength defect inspection
- Illumination emulates EUVL stepper
- Spatial resolution: 55 nm






*F. Brizuela et al, Optics Letters vol. 34, 271, (2009)*



# 13.2 nm laser-based microscope for defect inspection in EUV lithography masks






$\lambda = 13.2$  nm resonant with Mo/Si coatings in extreme ultraviolet lithography masks





EUV Optics from CXRO, Berkeley

## EUV MASK IMAGES

Spatial resolution: 55 nm half-pitch


 $\text{LER} < 10\% \text{ CD}$ *F. Brizuela et. al., Opt. Exp. 14, 14467 (2010)*



## Self-seeded

EUV Amplifier


Spontaneous emission



## Injection-seeded

EUV Amplifier

Coherent seed




Seed pulses can be greatly amplified preserving or even improving their properties




# Injection-seeding SXR Lasers have full phase-coherence and shorter pulsewidth



## Full spatial coherence




## Full temporal coherence



## Shorter pulsewidth ( $1.13 \pm 0.47$ )ps




# High repetition rate, diode-pumped table-top soft x-ray lasers



## Laser Diode Pumping Advantages





- Highly efficient
  - >50% Electrical efficiency
- Narrow bandwidth
  - Efficiently pump a single transition
- Directional
  - End-pumping
- Very high average power
  - Allow high repetition rate
- Compact







- Absorption bands at InGaAs wavelengths
- Small quantum defect (<10%)
- Long lifetime for high energy storage

# Beyond 100 Hz Repetition Rate: 500 Hz Compact High Energy, High Power Pump Laser








# Demonstration of bright 18.9 nm laser at 100 Hz repetition rate



$\lambda = 18.9 \text{ nm}$  Laser Average Power: **0.2 mW**

*B. Reagan, et al., Opt. Lett. 37 (2012).*

# $\lambda = 18.9 \text{ nm}$ : $10^5$ Consecutive EUV Laser Shots at 100 Hz repetition rate




*B. Reagan et al. Optics Express, 21, 28320 (2013)*

# Extending diode-pumped lasers to $\lambda = 10.9$ nm

## Ni-like ions


44.956 39 Y 88.906 71 Lu 174.97 lawrencium	47.867 40 Zr 91.224 72 Hf 178.49 rutherfordium	50.942 41 Nb 92.906 73 Ta 180.95 dubnium	51.996 42 Mo 95.94 74 W 183.84 seaborgium	54.938 43 Tc [98] 75 Re 186.21 bohrium	55.845 44 Ru 101.07 76 Os 190.23 hassium	58.933 45 Rh 102.91 77 Ir 192.22 meitnerium	63.546 46 Pd 106.42 78 Pt 195.08 unnilium	68.923 47 Ag 107.87 79 Au 196.97 ununnium	72.64 48 Cd 112.41 80 Hg 200.59 ununium	74.922 49 In 114.82 81 Tl 204.38 ununquadium	78.972 50 Sn 118.71 82 Pb 207.2 ununquadium	74.922 51 Sb 121.76 83 Bi 208.98 ununquadium	78.972 52 Te 127.60 84 Po [209] ununquadium	126.905 53 I 126.90 85 At [210] ununquadium	131.29 54 Xe 131.29 86 Rn [222] ununquadium
---	---	---	--	---	---	--	--	--	--	---	--	---	--	--	--





- 0.2 mW at  $\lambda = 18.9$  nm
- 0.1 mW at  $\lambda = 13.9$  nm



# Extension of gain-saturated table-top SXRL to sub-10 nm wavelengths using lanthanide ions




# Gain-saturated sub-10 nm table-top lasers











# Lasing in transitions down to 7.36 nm

Nickel-like lanthanide ions  $4d^1S_0 - 4p^1P_1$

57	58	59	60	61	62	63	64	65	66	67	68	69	70	71
La	Ce	Pr	Nd	Pm	Sm	Eu	Gd	Tb	Dy	Ho	Er	Tm	Yb	Lu
Lanthanum 138.905	Cerium 140.12	Praseodymium 140.908	Neodymium 144.24	Promethium 145	Samarium 150.36	Europium 151.964	Gadolinium 157.25	Terbium 158.92534	Dysprosium 162.505	Holmium 164.93032	Erbium 167.259	Thulium 168.93401	Ytterbium 173.04	Lutetium 174.967
[Xe]5d <sup>1</sup>	[Xe]4f <sup>1</sup>	[Xe]4f <sup>2</sup>	[Xe]4f <sup>3</sup>	[Xe]4f <sup>4</sup>	[Xe]4f <sup>5</sup>	[Xe]4f <sup>6</sup>	[Xe]4f <sup>7</sup>	[Xe]4f <sup>7</sup> 5d <sup>1</sup>	[Xe]4f <sup>9</sup>	[Xe]4f <sup>10</sup>	[Xe]4f <sup>11</sup>	[Xe]4f <sup>12</sup>	[Xe]4f <sup>13</sup>	[Xe]4f <sup>14</sup> 5d <sup>1</sup>
5.0760	5.5367	5.473	5.5293	5.562	5.6437	5.6704	5.7488	5.8038	5.8389	5.9211	6.1077	6.1843	6.2542	5.4258



Gain saturated




Wavelength (nm)

Wavelength (nm)





# Extension of gain-saturated table-top SXRL to sub-7 nm wavelengths using lanthanide ions



# Gain saturation in Ni-like Gd at 6.86 nm


57	58	59	60	61	62	63	64	65	66	67	68	69	70	71
La	Ce	Pr	Nd	Pm	Sm	Eu	Gd	Tb	Dy	Ho	Er	Tm	Yb	Lu
Lanthanum 138.9055 [Xe]5s <sup>2</sup> 5p <sup>6</sup>	Cerium 140.118 [Xe]4f <sup>1</sup> 5s <sup>2</sup> 5p <sup>6</sup>	Praseodymium 140.90765 [Xe]4f <sup>2</sup> 5s <sup>2</sup> 5p <sup>6</sup>	Neodymium 144.24 [Xe]4f <sup>3</sup> 5s <sup>2</sup> 5p <sup>6</sup>	Promethium (145) [Xe]4f <sup>4</sup> 5s <sup>2</sup> 5p <sup>6</sup>	Samarium 150.36 [Xe]4f <sup>6</sup> 5s <sup>2</sup> 5p <sup>6</sup>	Europium 151.96 [Xe]4f <sup>7</sup> 5s <sup>2</sup> 5p <sup>6</sup>	Gadolinium 157.25 [Xe]4f <sup>7</sup> 5s <sup>2</sup> 5p <sup>6</sup>	Terbium 158.92534 [Xe]4f <sup>9</sup> 5s <sup>2</sup> 5p <sup>6</sup>	Dysprosium 162.500 [Xe]4f <sup>10</sup> 5s <sup>2</sup> 5p <sup>6</sup>	Holmium 164.93032 [Xe]4f <sup>11</sup> 5s <sup>2</sup> 5p <sup>6</sup>	Erbium 167.259 [Xe]4f <sup>12</sup> 5s <sup>2</sup> 5p <sup>6</sup>	Thulium 168.93401 [Xe]4f <sup>13</sup> 5s <sup>2</sup> 5p <sup>6</sup>	Ytterbium 173.04 [Xe]4f <sup>14</sup> 5s <sup>2</sup> 5p <sup>6</sup>	Lutetium 174.967 [Xe]4f <sup>14</sup> 5s <sup>2</sup> 5p <sup>6</sup>



Energy on target = 14 J , pre-pulse width= 185 ps, short pulse width= 0.7 ps , 40% in pre-pulse




# Laser amplification down to 5.8 nm (Ni-like Dy)

57	58	59	60	61	62	63	64	65	66	67	68	69	70	71
La	Ce	Pr	Nd	Pm	Sm	Eu	Gd	Tb	Dy	Ho	Er	Tm	Yb	Lu
Lanthanum 138.9055 [Xe]5s <sup>2</sup> 5p <sup>6</sup> 6s <sup>2</sup>	Cerium 140.118 [Xe]4f <sup>1</sup> 5s <sup>2</sup> 5p <sup>6</sup>	Praseodymium 140.90765 [Xe]4f <sup>3</sup> 5s <sup>2</sup> 5p <sup>6</sup>	Neodymium 144.24 [Xe]4f <sup>4</sup> 5s <sup>2</sup> 5p <sup>6</sup>	Promethium (145) [Xe]4f <sup>5</sup> 5s <sup>2</sup> 5p <sup>6</sup>	Samarium 150.36 [Xe]4f <sup>6</sup> 5s <sup>2</sup> 5p <sup>6</sup>	Europium 151.964 [Xe]4f <sup>7</sup> 5s <sup>2</sup> 5p <sup>6</sup>	Gadolinium 157.25 [Xe]4f <sup>7</sup> 5s <sup>2</sup> 5p <sup>6</sup>	Terbium 158.925 [Xe]4f <sup>9</sup> 5s <sup>2</sup> 5p <sup>6</sup>	Dysprosium 162.50 [Xe]4f <sup>10</sup> 5s <sup>2</sup> 5p <sup>6</sup>	Holmium 164.93032 [Xe]4f <sup>11</sup> 5s <sup>2</sup> 5p <sup>6</sup>	Erbium 167.259 [Xe]4f <sup>12</sup> 5s <sup>2</sup> 5p <sup>6</sup>	Thulium 168.93421 [Xe]4f <sup>13</sup> 5s <sup>2</sup> 5p <sup>6</sup>	Ytterbium 173.04 [Xe]4f <sup>14</sup> 5s <sup>2</sup> 5p <sup>6</sup>	Lutetium 174.967 [Xe]4f <sup>14</sup> 5s <sup>2</sup> 5p <sup>6</sup>





# Acknowledgements





**Andrew Duffin (Pacific Northwest Laboratory)**

**Elliot Bernstein (CSU Chemistry Department)**

**Zone plates: Weilun Chao, Patrick Naulleau (Lawrence Berkeley National Laboratory)**

**Funding: National Science Foundation and US Department**


# Demonstration of gain saturation in the 7.3 nm line of Ni-like Sm




Energy on target = 11.4 J , pre-pulse width= 185 ps, short pulse width= 0.7 ps , 30% pre-pulse

# Frequency multiplication of periodic patterns by fractional Talbot lithography


Frequency multiplication factors up to 5X were demonstrated




(a)



(b)




(c)



(d)




# $\lambda = 18.9\text{nm}$ Laser operating at 400 Hz




400 Hz repetition rate

Raw EUV photodiode output of first 400 shots (1 second).



Single shot EUV spectrum showing 18.9nm Laser at 400 Hz repetition rate




C. Baumgarten, et al., Opt. Lett. **41** (2016).

# Recent research has shrunk capillary discharge SXRL to desk-top size



Smallest SXRL laser ,  $\lambda=46.9$  nm

12 Hz repetition rate, 0.15 mW average power



- 10 microjoule /pulse
- 0.15 mW average power
- 1-12 Hz repetition rate
- Pulse duration ~1.5 ns
- $\Delta\lambda/\lambda < 1 \times 10^{-4}$



S. Heinbuch, M. Grisham, D. Martz, J.J. Rocca  
*Optics Express*, 30,2095, (2005)






# Molecular sensitivity to aminoacid Alanine of 0.01 amol is 40× that of SIMS TOF

10 keV Bi<sup>+</sup> at 45° into C




ALANINE	SXR TOF (First results)	SIMS TOF
Area probed (μm <sup>2</sup> )	4.5 10 <sup>-2</sup>	2.25 10 <sup>4</sup>
Depth probed (μm)	3.5 10 <sup>-3</sup>	2 10 <sup>-3</sup>
Volume probed (aL)	0.05	4.5 10 <sup>4</sup>
Sensitivity (amol)	<b>0.01</b>	0.4
Ion yield, normalized to probed volume (L <sup>-1</sup> )	4 10 <sup>14</sup>	8.5 10 <sup>10</sup>
Level of fragmentation	1.1	1

# Electron impact excitation of 6.86 nm Gd laser requires plasma with high electron temperature






## Electron impact excitation rate $4d^1S_0$




## Previous work achieved unsaturated lasing at 6.86 nm in Ni-like Gd








•Daido et al. using 250 J of laser pump energy (Phys. Rev. Lett. 75 1074, 1995)



# Compact plasma-based soft x-ray lasers can be installed at the application's site



- High pulse energy ( $\mu\text{J}$ - $\text{mJ}$ )
- High monochromaticity ( $\lambda/\Delta\lambda < 10^{-4}$ )
- High peak spectral brightness











# Optimum grazing incidence angle increases linearly with Z

Measurement is indicative of the electron density where maximum gain occurs

Experimentally determined optimum angles up to Dy (Z=66)





# Plasma formation controlled pre-pulse duration, delay, is critical for large Sm gain



# Capillary plasma columns generate gain-saturated soft X-ray amplification in $\text{Ar}^{+8}$


Exponential amplification in  $\text{Ar}^{+8}$   $3p^1S_0$ - $3s^1P_1$  line at 46.9 nm




\*Phys. Rev. Lett. **73** 2192 (1994)

3-4 mm Capillary diameter, 25-40 kA, >30 ns rise time  
Electron Temperature: 60-100 eV, Electron Density:  $0.2\text{-}1.0 \times 10^{19} \text{ cm}^{-3}$

# UV PI on Uranium (~1%) filter glass shows higher SUE




26.4 eV SPI




- 6.19 eV – 1<sup>st</sup> ionization potential of U

UV PI


Post-ionization: 3.5 eV









<http://en.wikipedia.org/wiki/File:FEL.png>



M.Beyer et al. PNAC, (2010)





A. Barty et al. Nat. Phot., (2008)



Nagler et al. Nat. Phys. (2011)





# 0.014% sample utilization efficiency (SUE) of trace elements is similar to SIMS TOF



- 500 ppm of Uranium in NIST 611 glass
- 1200 shots → 10×10×0.3 μm crater


*T. Green et al, J. of Analytical Atomic Spectrometry, 32, 1092, (2017)*

## 3-D maps of materials composition with nanometer resolution






# Applications in dense plasma diagnostics and photochemistry




## Plasma Interferometry


J. Filevich et al PRL 94, 035005 (2005)



M. Purvis et al. Phys. Rev.E, **76**, (2007); **124**, (2010)



## Single photon ionization mass spectrometry




F. Dong et al. J.Chem.Phys **124**, (2006)

F. Dong et al. J.Am.Chem Soc. **131**, (2009)





# Scaling to shorter wavelengths requires hotter-denser plasmas




# Soft X-ray lasers excited by rapid heating of plasmas with short laser pulses

## Laser Pumping Geometry



Grazing incidence allows for efficient heating of plasma region with optimum electron density





R. Keenan et al, *Phys. Rev. Lett.* 94, 103901 (2005) ; B.M. Luther et al, *Opt. Lett.* 30, 165 (2005);  
Transient excitation: P. Nickels, V. Shlyaptsev et al. *Phys. Rev.Lett.* 78,2 748, (1997)

# Simulation showed gain-saturated amplification at 13.2 nm in Ni-like Cd can be achieved with ~ 1 J pump





Pre-pulse  
300 mJ, 120 ps

Heating pulse  
1 J, 6 ps







# Lasers pumped by a 5-10 Hz ~ 1 J Short Pulse Table-top Ti: Sapphire System



# High repetition rate table-top SXRL in transitions of Ni-like ions down to 10.9 nm





Gain saturated  
operation  
demonstrated



Y. Wang et al, *Phys. Rev. A* **72**, 053807 (2005)


\*D. Martz et al. *Optics Lett.* **35**, 1632 (2010)



## Self-seeded

EUV Amplifier


Spontaneous  
emission



## Injection-seeded

EUV Amplifier

Coherent  
seed




Seed pulses can be greatly amplified preserving or even improving their properties




# Injection-seeding SXR Lasers have full phase-coherence and shorter pulsewidth



## Full spatial coherence



## Full temporal coherence






## Shorter pulsewidth ( $1.13 \pm 0.47$ )ps



# 13.2 nm laser-based microscope for defect inspection in EUV lithography masks




$\lambda = 13.2$  nm resonant with Mo/Si coatings in extreme ultraviolet lithography masks




EUV Optics from CXRO, Berkeley


# Extension of gain-saturated table-top SXRL to sub-10 nm wavelengths using lanthanide ions



# Electron impact excitation of 8.8 nm La laser requires plasma with high electron temperature



## Electron impact excitation rate $4d\ 1S_0$





## Previous work achieved unsaturated lasing at 8.8 nm in Ni-like La

• *Daido et al. using 520 J of laser pump energy (Optics Lett. 21, 958, 1996)*





*Kawachi et al. using 18 J picosecond pulses (Phys. Rev. A, 69, 2004)*

Gekko XII Laser (Osaka)









# Simulation for 8.8 nm table-top Laser in Ni-like La predicts $< 7$ J pump energy needed for gain saturation





Simulation by Mark Berrill

# Titanium-Sapphire pump laser


Average Energy Pre-compression= 13 J  
Std div. = 1.5 %








# Gain-saturated sub-10 nm table-top lasers






# Demonstration of Gain-saturated table-top laser at 8.8 nm at 1 Hz repetition rate

Ni-like Lanthanum  $4d^1S_0 - 4p^1P_1$




7.5 J Total Pump Energy

Pulse energy up to  $\sim 2.7 \mu\text{J}$





# 1 Hz $\lambda = 8.8$ nm laser output intensity exceeds computed saturation intensity by an order of magnitude

## Near field beam profile measurement



## 1 Hz repetition rate




SXRL Fluence:  $0.6 \text{ J cm}^{-2}$   
 Experiment:  $I \sim 2.4 \times 10^{11} \text{ W cm}^{-2}$   
 Computed  $I_{\text{sat}}: \sim 3 \times 10^{10} \text{ W cm}^{-2}$


*D. Alessi et al. Phys. Rev. X, 1, 021023 (2011)*

# Lasing in transitions down to 7.36 nm

## Nickel-like lanthanide ions $4d^1S_0 - 4p^1P_1$



# Gain-saturated table-top SXRLs cover 8.8 nm - 47 nm wavelength region






# The Next Generation: Increasing the repetition rate of Table-Top Soft X-Ray Lasers to 100 Hz


Laser Diode Drivers



Solid State Ultrashort Pulse High Power Laser




Soft X-Ray Plasma Amplifier



13.9 nm  
laser

$\text{Ag}^{+19}$

e




# Directly diode-pumped Yb CPA laser increases repetition rate and average power

## Laser Diode Pumping Advantages




- Highly efficient
  - >50% Electrical efficiency
- Narrow bandwidth
  - Efficiently pump a single transition
- Directional
  - End-pumping
- Very high average power
  - Allow high repetition rate
- Compact



- Absorption bands at InGaAs wavelengths
- Very low quantum defect (<10%)
- Long lifetime for high energy storage

# Thermal and gain properties of Yb:YAG are dramatically improved at cryogenic temperature

Yb:YAG at room and cryogenic temperature	300 K	77 K	
Thermal conductivity (W/mK)	10	90	x9
Thermo-optic coefficient ( $10^{-6}/K$ )	7.8	0.9	x1/7
Expansion coefficient ( $10^{-6}/K$ )	6.14	1.95	x1/4
Saturation fluence ( $J/cm^2$ )	9.2	1.7	x1/7




G. A Slack and D. W. Oliver; *Phys. Rev. B* **4**; 592-609 (1971)  
 R. Wynne, J. L. Daneu and T. Y. Fan; *Appl. Opt.* **38**, 3282-3284 (1999)  
 R.L. Aggarwal, et. al., *Journal of Applied Physics*, **98**, 103514, (2005).

## Other recent cryogenic diode-pumped CPA work:

1. K.H. Hong, et al., *Optics Letters* **35**, 1752, (2010).
2. D. Rand, et al., CM3D.4 CLEO 2012.
3. D.E. Miller, et al., CM3D.2 CLEO 2012.
4. K. Ogawa, et al., CMB.4, CLEO 2011.



# Compact high power diode-pumped CPA laser driver for 100 Hz table-top SXRL








# 2<sup>nd</sup> stage cryo-cooled Yb:YAG amplifier

140 mJ, 100 Hz, amplifier operation demonstrated


## 2nd Stage Amplifier




Single pass gain





Pulse energy




Beam pattern







*Uncompressed pulses*




*$M^2$  of amplified pulses*




*2<sup>nd</sup> order autocorrelation of compressed 1 J pulses*



**1 J, 5 ps pulses at 100 Hz repetition rate**




# 100 Hz Operation






# Gain-Saturated 18.9nm Laser Operation at 100 Hz repetition rate



940 mJ on target target moved at 200  $\mu\text{m/s}$ , (2 $\mu\text{m}/\text{shot}$ )




**Mean Energy = 1.46  $\mu\text{J}$ ,  $\sigma = 11\%$**


**0.15 mW average power**

( Fermi FEL 20-65 nm: 30-60  $\mu\text{J}$  x 10 Hz = 0.3-0.6 mW Luca Giannessi ICXRL)

# Helicoidal targets developed to allow continuous operation at 100 Hz repetition rate



*Slab targets for parameterization of the soft x-ray laser*




Infrared Laser  
Pulses

Soft X-Ray  
Laser

*Helicoidal target for applications*


# Demonstration of *all-diode-pumped laser* at 13.9nm in Ni-like silver plasma





# Summary

- Gain-saturated table-top SXRLs reach  $\lambda = 8.85$  nm. Amplification observed down to  $\lambda = 7.3$
- Compact diode-pumped soft x-ray laser operating at record 100 Hz rep. rate produces 0.15 mW average power on a table-top



Work Supported by the NSF Engineering Research Centers Program  
and the US Department of Energy



# Acknowledgement





# New Facilities are allowing to expand optical drivers for secondary radiation generation

Advanced  
Beam  
Laboratory



Diode Pumped Ti:Sa CPA



Diode Pumped Yb:YAG CPA

# Photonics at extreme ultraviolet and soft x-ray wavelengths on a table-top

**Carmen S. Menoni**

*Center for Extreme Ultraviolet Science and Technology*

*Department of Electrical & Computer Engineering*


*Department of Chemistry*

**IEEE Distinguished Lecture**


**Colorado  
State  
University**




**Imaging of nanostructures**



**Movies at the nanoscale**




**Molecular Imaging**





# EUV/SXR sources are paving the way to scientific research and technological innovation

## PHOTON FACTORIES




## COMPACT SOURCES


### Incoherent sources







### High Harmonic Generation






### Lasers



# Compact plasma-based soft x-ray lasers can be installed at the application's site



- High pulse energy ( $\mu\text{J}$ -mJ)
- High monochromaticity ( $\lambda/\Delta\lambda < 10^{-4}$ )
- High peak spectral brightness





# New Facilities are allowing to expand optical drivers for secondary radiation generation

Advanced  
Beam  
Laboratory




Diode Pumped Ti:Sa CPA

Diode Pumped Yb:YAG CPA



# Gain-saturated table-top SXRLs cover 8.8 nm - 47 nm wavelength region - Pump: CPA Ti:Sapphire






# SXRL microscopes are critical photonic technologies for imaging nanostructures and surfaces

## Transmission


$\lambda=46.9$  nm, 13.9 nm



38 nm spatial resolution @ 13.9 nm

*G. Vaschenko et al, Opt.Lett 2006*

## Single shot Imaging of nanostructures




50 nm Carbon Nanotubes

Spatial resolution: 50 nm

*C. Brewer et al, Opt.Lett 2008*

## Time resolved Imaging




Frames from a SXR motion picture


*S. Carbajo et al, Opt. Lett . 2012*

## Reflection

$\lambda=46.9$  nm, 13.2 nm



## Aerial inspection of EUVL masks ( $\lambda=13.2$ nm)




*F. Brizuela et al, Op.Lett. 2009, Opt. Exp. 2010*

# Compact $\lambda = 46.9$ nm full field microscope

## Laser Output Characteristics:


- Output Energy:  $> 10 \mu\text{J}/\text{pulse}$
- Repetition rate: up to 12 Hz
- Pulse duration:  $\sim 1.5$  ns
- Bandwidth:  $\Delta\lambda/\lambda < 1 \times 10^{-4}$
- Tailored spatial coherence

Sc/Si multilayer coated  
Schwarzschild condenser  
Throughput: 13%  
NA=0.18,  
working distance: 60 mm




Freestanding Fresnel  
zone plate objective

NA	$\Delta r$	f
0.32	70 nm	0.75 mm
0.19	120 nm	1.28 mm



Capillary discharge pumped  
 $\lambda = 46.9$  nm laser



vacuum chamber

EUV- sensitive CCD  
connected to computer

Power supply and  
gas handling controls


Courtney  
Brewer

# EUV/SXR microscope captures images of nanostructures with very high resolution

## TRANSMISSION


NA=0.32 (M~1000)

200 nm half period diatom



5 sec exposure

50 nm carbon nanotube




Single shot,  
1.5 ns  
exposure

C. Brewer, et al, Opt. Lett. 33, 518 (2008)


## REFLECTION

NA=0.12 and 0.19 (M~250)


5-20 sec exposures



Partially Processed  
semiconductor chip



100 nm thick GaN nanowire  
between Al contacts



Zr surface showing twin in grain

F. Brizuela et al, Opt.  
Express 13 :3983, (2005).



# Single-shot sequential imaging



## Motion picture of nanoscale dynamics

**Laser Output at 46.9 nm**

**Energy > 10  $\mu\text{J}/\text{pulse}$**   
**( $2.4 \cdot 10^{12}$  ph/pulse)**


**Pulse duration ~1 ns**

- Test object: Oscillating cantilever driven by periodic voltage signal
- Single shot images were acquired to map oscillation over an entire period







# Visualizing magnetic nanoscale dynamics



Experiments






Simulations



- Variations of 30 nm in the amplitude of the oscillation detected.
- Model predicts change in the restoring force of the tip and associated amplitude changes that agree with experiment

# SXRL laser holography captures images with a single laser flash



## IN LINE HOLOGRAPHY



**50 nm spatial resolution**

*P. Wachulak, JOSA B, 25, 1811(2008).*

## SINGLE SHOT EUV FOURIER HOLOGRAPHY CAPTURES MOTION OF NANOPILLARS




**80 nm spatial resolution**

*E. Malm et al, Express:21:9959 (2013)*

**Marconi's group**


# EUV/SXR interference contrast imaging




J. Nedjl et al, submitted to IEEE Photonics Journal

# Full field images of a 230 nm dense line elbow pattern with 600x magnification

## RAW IMAGE




## PHASE IMAGE




Phase difference of  $(2.3 \pm 0.3)$  rad corresponds to  $(100 \pm 15)$  nm thick Si

## ABSORPTION IMAGE





## PHASE CUTS





# $\lambda=13.2$ nm microscope for EUVL mask inspection

## illumination emulates EUVL stepper




# Nanoscale Molecular imaging by SXR laser ablation mass spectrometry


**Mass Spectrometry Imaging (MSI) allows for identification of and mapping of the spatial distribution of elemental and molecular components in solid samples**

CHALLENGES IN MSI of bio/organic samples to reach nanoscale resolution are:

- Reduction of ion yield when probing nanoscale volumes
- Molecular fragmentation




# Concept of EUV/SXR Mass Spectrometry Imaging



- Laser pulses are focused onto specimen
  - Wavelength: 46.9 nm
  - Pulse duration: 1.5 ns
  - Pulse energy: 0.01 mJ
- The focused laser ablates the material and simultaneously ionizes the fragments in the ablation plume
- Ions are extracted into the Time of Flight mass spectrometer

# $\lambda=46.9$ nm laser ablation offers unique characteristics for MSI


- **Exceptional lateral resolution:**  $\lambda = 46.9$  nm light can be focused to  $\sim 100$  nm spots [G. Vaschenko et al, Opt. Lett. Vol. 31, 3615-3617]
- **Exceptional depth resolution:**  $\lambda = 46.9$  nm light is strongly absorbed in most materials
- **EUV photons break bonds in organic materials** [L. Juha et al, APL, Vol. 89, 034109, 3009]



Profile of a crater with a volume of 1 aL ablated in PMMA with a single shot. The crater's profile shows no sign of thermal damage


## EUV laser ablation and plasma generation

- Plasma is transparent to EUV photons
- Electron heating by inverse bremsstrahlung is negligible
- EUV photons absorbed by photoionization have sufficient energy to single photon-ionize any atom or molecule






# 3D imaging EUV mass spectrometry imaging nanoprobe



**EUV Optics**  
Fresnel zone plate 0.16 NA



## **EUV laser**

- Wavelength: 46.9 nm
- Energy per pulse > 10 μJ
- Repetition rate: 12 Hz
- Pulse duration ~1.5 ns

*S. Heinbuch Optics Express  
vol. 13, 4050 (2005)*


Optics engineered by  
W. Chao and E. Anderson at  
Center of X-Ray Optics,  
Lawrence Berkeley Lab.

*E.H. Anderson, IEEE J. Quantum  
Electronics, vol.42, 27, 2006*

## Results:




# Molecular imaging of inorganic samples (negative mode extraction)

Mass spectrum of SiO<sub>2</sub>



# Results:


## Depth profiling in organic samples with resolution $\sim 20$ nm




## Results:

# 2D imaging with 70 nm lateral resolution


2D ion image plotting content of resist determined from the intensity average of 4 peaks in the 70-120 m/z range



AFM image of ablated region



Edge profile from AFM







# Results:

## 3D composition mapping of a single *M. Smegmatis* bacterium

### Iso-lines of two significant lipid fragments






**Voxel size:  $0.3 \times 0.3 \times 0.08 \mu\text{m}^3$**

*I. Kuznetsov, et al, to be published  
in Nat. Comm. 2015*

# Comparison of EUV TOF with leading molecular mass spectrometry methods UV LDI TOF and SIMS TOF

## Mass spectra of alanine




- Strong absorption
- High 3D localization
- Efficient photo-ionization
- **Sensitivity: 0.01 amol**
- Fragmentation: 1.1

- Negligible absorption
- Absence of ionization
- (Requires Matrix)

- Ion collision process
- **Sensitivity: 0.4 amol**
- Fragmentation: 1

# SXR Coherent lithography prints arbitrary motifs defect free

Talbot Distance




- NON- CONTACT
- LENSLESS
- LITHOGRAPHY
- SHORT EXPOSURES, 20-50 LASER SHOTS
- PRINTS PATTERNS DEFECT FREE OVER AREAS  $0.5 \times 0.5 \text{ mm}^2$



*L. Urbanski, Opt. Lett. 37, 3633 (2012).*

# Error-Free Printing of Periodic metallic structures

Talbot Mask



Print in HSQ




**Grating etched in Au  
500/600 lines/spaces**



# Smaller Features: Fractional Talbot EUV Lithography

Periodic masks produce images with higher frequency in positions that are a fraction of the Talbot distance




Talbot carpet

Intensity cross sections


Optics Letters, Doc. ID 224904

# Frequency multiplication of periodic patterns by fractional Talbot lithography


Frequency multiplication factors up to 5X were demonstrated




(a)



(b)



(c)



(d)

## Fractional Talbot lithography with extreme ultraviolet light


Hyun-su Kim<sup>1,4,\*</sup>, Wei Li<sup>3</sup>, Serhiy Danylyuk<sup>2</sup>, William S. Brocklesby<sup>4</sup>,  
Mario C. Marconi<sup>3</sup>, and Larissa Juschkina<sup>1</sup>

AFM measurements for the lithography results, showing (a)  $M_{sf} = 1$ , (b)  $M_{sf} = 2$ , (c)  $M_{sf} = 3$ , and (d)  $M_{sf} = 5$  from the parent mask of 3  $\mu\text{m}$  pitch.

Overlaid red curves show corresponding AFM cross-sections. The height of the PMMA structures (a-d) was  $\sim 80$  nm.


*Optics Letters, Doc. ID 224904*

# Double exposure and stitching prints 40 nm lines over millimeter square areas




**Smallest lines printed**

**40 nm / 110 nm lines / Spaces**



**Coherent Printing with 100 Hz,  
18.9 nm laser**

**400 nm lines**



*1 Minute Exposure,  
0.5 mm<sup>2</sup> Patterned  
Area*

**Potential to print  
10 nm lines**

Non contact printing method is  
well suited for studies of EUV  
resist response

L. Urbanski et al

# Thanks to Collaborators

## Microscopy Team

J. Nedjl, S. Carbajo, I. Howlett, N. Monserut, E. Malm, M.C. Marconi and K. Buchanan (CSU) (recent work)

*Pioneers: G. Vaschenko, F. Brizuela*

## Mass spectrometry imaging team

I. Kuznetsov, J. Filevich, M. Woolston, T. Green, E. R. Bernstein, D.C. Crick, and J.J. Rocca (CSU)

## Nanopatterning Team

W. Li, L. Urbanski, M.C. Marconi

## Laser Engineering Team

Y. Wang, B. Reagan, M. Woolston, C. Baumgartner, A. Rockwood, S. Wang, L. Yin, J. Rocca (CSU)

## Optics Engineering Teams

W. Chao, E. H. Anderson, K. Goldberg, P. Naulleau

*Center for X-Ray Optics, Lawrence Berkeley Lab*

A. Vinogradov, I. Artioukov

*Levedev Physical Institute, Russia*




**AIR Research Alliance project “  
*Development of Key Technology for next generation  
projection lithography of integrated circuits at 6.X  
nm wavelengths”***

**Collaboration between EUV ERC and Cymer**


**Funded by NSF and Cymer**

# BRIGHT SXR LASERS ENABLE NOVEL TOOLS FOR NANOSCIENCE AND NANOTECHNOLOGY ON A TABLE TOP


## Imaging of nanostructures




## Movies of nanoscale interactions




## 3D Molecular Imaging



## MICROSCOPY



## 3D MOLECULAR IMAGING



FUTURE

TODAY





**Rocky Mountain National Park  
FALL**



**Sky Colorado - WINTER**



INTERNATIONAL  
YEAR OF LIGHT  
2015

**Thank you!**

conscious community



**Horsetooth Park , Fort Collins -  
SUMMER**